

The All New iZeePulse Home Unit

THE PERFECT PEMF HOME UNIT

icbfamily of companies



Solid State digital iZeePulse Home/Clinic Unit

Length: 10.25 in, Width: 12.5 in, Height: 4 in, Weight: 5.8 lbs VOLTAGE:120 to 230 VAC, 50 to 60H Gauss:7,500G | TESLA: .75T

Pulsed Electro Magnetic Field Therapy



Length: 7.5 Inch, Width: 10.5 Inch, Height: 3 Inch, Weight: 4.4 lbs VOLTAGE:120 to 230 VAC, 50 to 60Hz Gauss:7,500G | TESLA: .75T







\$8,495 ATTACHMENTS: \$8,495 DOUBLE LOOP (7" DIAMETER) Suggested Retail Price MAT (18" X24") -\$500 (Instant Discount) \$7,995

Special Price

YOUR CHOICE: DESKTOP MINI TOTE

INFO@IZEEPULSE.COM

tebramilyorcompanies

Difference between Hummers & Ringers

"Hummers"

These units use a low intensity gauss PEMF to deliver preset frequencies that **resonate receptors** at the cell membrane

PEMF vibration is NOT felt by the person

Requires longer treatment time sessions

"Ringers" Our PEMF Devices

These units use a high intensity PEMF Pulse to deliver harmonic impulses that **resonate the entire cell membrane**.

Instant results and mild to strong PEMF sensations are felt by the person

Short term treatment time sessions can achieve results

INFO@IZEEPULSE.COM




A sample of results achieved with the home/Clinic unit

My Surgeon was amazed on how fast I healed!!! This healing happened in less than 28 days thanks to the iZeeFulse Home **SEMF Unit**

ic bfamily of companies







\$7,995 iZeePulse Home/Clinic Unit

GET BACK WITH THE PERSON WHO REFERRED YOU TO THIS PRESENTATION TO GET YOUR QUESTIONS ANSWERED OR YOU CAN CALL CORPORATE (563) 726-0120

CebiamIlyacompanies

